

### **Multi-drug Rapid Test (Urine) Package Insert**

Version: Z

Effective Date: 2023-10-26

For professional in vitro diagnostic use only.

instruction Sheet for testing of any combination of the following drugs: AMP/BAR/BZO/BUP/COC/THC/MTD/MET/MDMA/MOP/OPI/PCP/TCA/TRA/KET/OXY/EDDP/FYL/

AMPIBARIBZO/BUPICOCTHE/MI ID/MELIMIDED STATE AND ADDRESS AND ADDRE

The Multi-Drug Rapid Test Panel is a rapid chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in urine at the following cut-off concentrations:

Test	Calibrator	Cut-off (ng/mL)
Amphetamine (AMP 300)	d-Amphetamine	300
Barbiturates (BAR 300)	Secobarbital	300
Benzodiazepines (BZO 300)	Oxazepam	300
Buprenorphine (BUP 10)	Buprenorphine	10
Cocaine (COC 150)	Benzoylecgonine	150
Marijuana (THC 50)	11-nor-Δ <sup>9</sup> -THC-9 COOH	50
Methadone (MTD 300)	Methadone	300
Methamphetamine (MET 300)	d-Methamphetamine	300
Methylenedioxymethamphetamine (MDMA 300)	d,l-Methylenedioxymethamphetamine	300
Morphine (MOP 300)	Morphine	300
Opiate (OPI 2000)	Morphine	2000
Phencyclidine (PCP 25)	Phencyclidine	25
Tricyclic Antidepressants (TCA 1,000)	Nortriptyline	1,000
Tramadol (TRA 100)	Cis-Tramadol	100
Ketamine (KET 300)	Ketamine	300
Oxycodone (OXY 100)	Oxycodone	100
2-ethylidene-1,5-dimethyl-	2-ethylidene-1,5-dimethyl-	100
3,3-diphenylpyrrolidine (EDDP 100)	3,3-diphenylpyrrolidine	
Fentanyl(FYL200)	Fentanyl	200
Propoxyphene(PPX300)	Propoxyphene	300
Synthetic Cannabis(K2 50)	JWH-018 5-Pentanoic acid	50
	JWH-073 4-butanoic acid	50
Pregabalin(PGB 500)	Pregabalin	500
6-Monoacetylmorphine(6-MAM 10)	6-Monoacetylmorphine	10
Methcathinone(MCT 500)	S(-)-Methcathinone	500
Ethyl- β-D-Glucuronide(ETG 500)	Ethyl- β-D-Glucuronide	500

Configurations of the Multi-Drug Rapid Test Panel come with any combination of the above listed drug analytes. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when reliminary positive results are indicated

The Multi-Drug Rapid Test Panel is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to selectively detect elevated levels of specific drugs in urine. Amphetamine (AMP)

Amphetamine is a Schedule II controlled substance available by prescription (Dexedrine®) and is also available on the illicit market. Amphetamines are a class of potent sympathomimetic agents with therapeutic applications. They are chemically related to the human body's natural catecholamines; epinephrine and norepinephrine. Acute higher doses lead to enhanced stimulation of the central nervous system (CNS) and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. Cardiovascular responses to amphetamines include increased blood pressure and cardiac arrhythmias. More acute responses produce anxiety, paranoia, hallucinations, and psychotic behavior. The effects of Amphetamines generally last 2-4 hours following use and the drug has a half-life of 4-24 hours in the body. About 30% of eted in the urine in unchanged form, with the remainder as hydroxylated and deaminated deri The Multi-Drug Rapid Test Panel yields a positive result when the concentration of amphetamines in urine exceeds

### Barbiturates (BAR)

Barbiturates are CNS depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of barbiturates leads to tolerance and physical dependence.

Short-acting barbiturates taken at 400 mg/day for 2-3 months can produce a clinically significant degree of physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death.

Only a small amount (less than 5%) of most barbiturates are excreted unaltered in the urine. The approximate detection time limits for barbiturates are

Short acting (e.g. Secobarbital)	100 mg PO (oral)	4.5 days
Long acting (e.g. Phenobarbital)	400 mg PO (oral)	7 days <sup>2</sup>

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of barbiturates in urine exceeds detective

### Benzodiazepines (BZO)

zodiazepines are medications that are frequently prescribed for the symptomatic treatment of anxiety and sleep

disorders. They produce their effects via specific receptors involving a neurochemical called gamma aminobutyric acid (GABA). Because they are safer and more effective, benzodiazepines have replaced barbiturates in the treatment of both anxiety and insomnia. Benzodiazepines are also used as sedatives before some surgical and medical procedures, and for the treatment of seizure disorders and alcohol withdrawal.

treatment of seizure disorders and alcohol withdrawal.

Risk of physical dependence increases if benzodiazepines are taken regularly (e.g., daily) for more than a few months, especially a higher than normal doses. Stopping abruptly can bring on such symptoms as trouble sleeping, gastrointestinal upset, feeling unwell, loss of appetite, sweating, trembling, weakness, anxiety and changes in perception.

Only trace amounts (less than 1%) of most benzodiazepines are excreted unaltered in the urine; most of the concentration in urine is conjugated drug. The detection period for benzodiazepines in urine is 3-7 days.

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of benzodiazepines in urine exceeds

### Buprenorphine (BUP)

Buprenorphine (BUP)

Buprenorphine is a potent analgesic often used in the treatment of opioid addiction. The drug is sold under the trade names Subutex™, Buprenex™, Temgesic™ and Suboxone™, which contain Buprenorphine HCl alone or in combination with Naloxone HCl. Therapeutically, Buprenorphine is used as a substitution treatment for opioid addicts. Substitution treatment is a form of medical care offered to opiate addicts (primarily heroin addicts) based on a similar or identical substance to the drug normally used. In substitution therapy, Buprenorphine is as effective as Methadone but demonstrates a lower level of physical dependence. Concentrations of free Buprenorphine and Norbuprenorphine in urine may be less than 1 ng/ml after therapeutic administration, but can range up to 20 ng/ml in abuse situations. The plasma half-life of Buprenorphine is 2-4 hours. While complete elimination of a single dose of the drug can take as long as 6 days, the window of detection for the

parent drug in urine is thought to be approximately 3 days.

Substantial abuse of Buprenorphine has also been reported in many countries where various forms of the drug are available. The drug has been diverted from legitimate channels through theft, doctor shopping, and fraudulent prescriptions, and been

abused via intravenous, sublingual, intranasal and inhalation routes.

The Multi-Drug Rapid Test Panel yields a positive result when the Buprenorphine in urine exceeds detective level.

Cocaine is a potent central nervous system stimulant and a local anesthetic. Initially, it brings about extreme energy and

Cocaine is a potent central nervous system is introduced in a discontruction and a cocain and a cocain accusate, in many, it orings about extreme energy are restlessness while gradually resulting in tremors, over-sensitivity and spasms. In large amounts, occaine causes fever, unresponsiveness, difficulty in breathing and unconsciousness.

Cocaine is often self-administered by nasal inhalation, intravenous injection and free-base smoking. It is exerted in the urine in a short time primarily as benzoylecgonine. "Benzoylecgonine, a major metabolite of cocaine, has a longer biological half-life (5-8 hours) than cocaine (0.5-1.5 hours), and can generally be detected for 24-48 hours after cocaine

exposure. The Multi-Drug Rapid Test Panel yields a positive result when the concentration of benzoylecgonine in urine exceeds Marijuana (THC)

THC (\( \frac{A}\) clerahydrocannabinol) is the primary active ingredient in cannabis (marijuana). When smoked or orally administered, THC produces cuphoric effects. Users have impaired short-term memory and slowed learning. They may also experience transient episodes of confusion and anxiety. Long-term, relatively heavy use may be associated with behavioral disorders. The peak effect of marijuana administered by smoking occurs in 20-30 minutes and the duration is 90-120 minutes after one cigarette. Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 3-10 days after smoking. The main metabolite excreted in the urine is 11-nor-Δ9-tetrahydrocannabinol-9-carboxylic acid

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of THC-COOH in urine exceeds detective

Methadone is a narcotic analgesic prescribed for the management of moderate to severe pain and for the treatment of opiato dependence (heroin, Vicodin, Percocet, morphine). The pharmacology of oral methadone is very different from IV methadone. Oral methadone is partially stored in the liver for later use. IV methadone acts more like heroin. In most states

you must go to a pain clinic or a methadone maintenance clinic to be prescribed methadone. Methadone is a long acting pain reliever producing effects that last from twelve to forty-eight hours. Ideally, methadone frees the client from the pressures of obtaining illegal heroin, from the dangers of injection, and from the emotional roller coaster that most opiates produce. Methadone, if taken for long periods and at large doses, can lead to a very long withdrawal priod. The withdrawals from methadone are more prolonged and troublesome than those provoked by heroin cessation, yet the substitution and phased removal of methadone is an acceptable method of detoxification for patients and

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of methadone in urine exceeds detec

Methamphetamine is an addictive stimulant drug that strongly activates certain systems in the brain. Methamphetamine is closely related chemically to Amphetamine, but the central nervous system effects of Methamphetamine are greater Methamphetamine is made in illegal laboratories and has a high potential for abuse and dependence. The drug can be taken orally, injected, or inhaled. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. Cardiovascular responses to Methamphetamine include increased blood pressure and cardiac arrhythmias. More acute responses produce anxiety paranoia, hallucinations, psychotic behavior, and eventually, depression and exhaustion

The effects of Methamphetamine generally last 2-4 hours and the drug have a half-life of 9-24 hours in the body Methamphetamine is excreted in the urine primarily as Amphetamine, and oxidized and deaminated derivatives. However 10-20% of Methamphetamine is excreted unchanged. Thus, the presence of the parent compound in the urine indicates Methamphetamine use. Methamphetamine is generally detectable in the urine for 3-5 days, depending on urine pH level. The Multi-Drug Rapid Test Panel is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Methamphetamine in urine. The Multi-Drug Rapid Test Panel yields a positive result when the Methamphetamine in urine exceeds detective level.

Methylenedioxymethamphetamine (MDMA) Methylenedioxymethamphetamine (ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the treatment of obesity.5 Those who take the drug frequently report adverse effects, such as increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is thought to be via release of the neurotransmitter serotonin. MDMA may also release dopamine, although the general opinion is that this is a secondary effect of the drug (Nichols and Oberlender, 1990). The most pervasive effect of MDMA, occurring in virtually all people

who took a reasonable dose of the drug, was to produce a clenching of the jaws.

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of Methylenedioxymethamphetamine in

### Morphine (MOP)

e refers to any drug that is derived from the opium poppy, including the natural products, morphine and codeine, and the semi-synthetic drugs such as heroin. Opioid is more general, referring to any drug that acts on the opioid receptor. Opioid analgesics comprise a large group of substances which control pain by depressing the CNS. Large doses of morphine can produce higher tolerance levels, physiological dependency in users, and may lead to substance abuse.

Morphine is excreted unmetabolized, and is also the major metabolic product of codeine and heroin. Morphine is detectable in the urine for several days after an opiate dose.

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of morphine in urine exceeds detective

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of morphine in urine exceeds 2000 ng/mL This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).1 See morphine (MOP 300) for summary.

Phencyclidine, also known as PCP or Angel Dust, is a hallucinogen that was first marketed as a surgical anesthetic in the 1950's. It was removed from the market because patients receiving it became delirious and experienced hallucinations or PCP is used in powder, capsule, and tablet form. The powder is either snorted or smoked after mixing it with marijuans or PCP either the properties of the propert orally. After low doses, the user thinks and acts swiftly and experiences mood swings from euphoria to depression.

Self-injurious behavior is one of the devastating effects of PCP.

PCP can be found in urine within 4 to 6 hours after use and will remain in urine for 7 to 14 days, depending on factors such as metabolic rate, user's age, weight, activity, and diet.6 PCP is excreted in the urine as an unchanged drug (4% to 19%) and conjugated metabolites (25% to 30%).6

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of phencyclidine in urine exceeds 25

ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA). Tricyclic Antidepressants (TCA)

# TCA (Tricyclic Antidepressants) are commonly used for the treatment of depressive disorders. TCA overdoses can result in The Articepressants are commonly used for the treatment of depressave disorders. It A overdoses can result profound CNS depression, cardiotoxicity and anticholinergic effects. TCA overdose is the most common case of death from prescription drugs. TCAs are taken orally or sometimes by injection. TCAs are metabolized in the liver. Both TCAs and their metabolities are excreted in urine mostly in the form of metabolities for up to ten days. The Multi-Drug Rapid Test Panel yields a positive result when the concentration of trievelic antidepressants in urine exceeds 1,000 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not

have a recommended screening cut-off for tricyclic antidepressant positive specimens.

Tramadol(TRA) is a quasi-narcotic analgesic used in the treatment of moderate to severe pain. It is a synthetic analog of codeine, but has a low binding affinity to the mu-opioid receptors. Large doses of tramadol can develop tolerance and physiological dependency and lead to its abuse. Tramadol is extensively metabolized after oral administration. Approximately 30% of the dose is exerted in the urine as unchanged drug, whereas 60% is excreted as metabolites. The major pathways appear to be N- and O- demethylation, glucoronidation or sulfation in the liver.

The Multi-Drug Rapid Test Panel is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Tramadol in urine. The Multi-Drug Rapid Test Panel yields a positive result when Tramadol in urine exceed detective level.

Ketamine is a dissociative anesthetic developed in 1963 to replace PCP (Phencyclidine). While Ketamine is still used in human anesthesia and veterinary medicine, it is becoming increasingly abused as a street drug. Ketamine is molecularly similar to PCP and thus creates similar effects including numbness, loss of coordination, sense of invulnerability, muscle rigidity, aggressive / violent behavior, slurred or blocked speech, exaggerated sense of strength, and a blank stare. There is depression of respiratory function but not of the central nervous system, and cardiovascular function is maintained. The effects of Ketamine generally last 4-6 hours following use. Ketamine is excreted in the urine as unchanged drug (2.3%) and

The Multi-Drug Rapid Test Panel is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Ketamine in urine. The Multi-Drug Rapid Tes Panel yields a positive result when Ketamine in urine exceeds detective level.

Oxycodone is a semi-synthetic opioid with a structural similarity to codeine. The drug is manufactured by modifying Oxycodone is a surrey-minete opine with a state and a similarity of executive. It is up is maintained by incarrying thebaine, an alkaloid found in the opium poppy. Oxycodone, like all opiate agonists, provides pain relief by acting on opioid receptors in the spinal cord, brain, and possibly directly in the affected tissues. Oxycodone is prescribed for the relief of moderate to high pain under the well-known pharmaceutical trade names of OxyContin®, Tylox®, Percodan® and

Percocet®. While Tylox®, Percodan® and Percocet® contain only small doses of oxycodone hydrochloride combined with other analgesics such as acetaminophen or aspirin, OxyContin consists solely of oxycodone hydrochloride in a time-release form. Oxycodone is known to metabolize by demethylation into oxymorphone and noroxycodone. In a 24-hour urine, 33-61% of a single, 5 mg oral dose is excreted with the primary constituents being unchanged drug (13-19%), conjugated drug (7-29%) and conjugated oxymorphone (13-14%). The window of detection for Oxycodone in urine is expected to be

similar to that of other opioids such as morphine.

The Multi-Drug Rapid Test Panel is a rapid urine screening test that can be performed without the use of an instrument. The

Interstulling a monoclonal antibody to selectively detect elevated levels of Oxycodone in urine. The Multi-Drug Rapid Test Panel syleids a positive result when Oxycodone in urine exceeds 100ng/ml.

2-ethyliden-1,5-dlmethyls-3,a-diphenylpyrrolldine (EDDP)

Methadone is an unusual drug in that its primary urinary metabolites (EDDP and EMDP) are cyclic in structure, making them very difficult to detect using immunoasays targeted to the native compound. "Exacerbating this problem, there is a subsection of the population classified as "extensive metabolizers" of methadone. In these individuals, a urine specimen may not contain enough parent methadone to yield a positive drug screen even if the individual is in compliance with their methadone maintenance. EDDP represents a better urine marker for methadone maintenance than unmetabolized

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of EDDP in urine exceeds detective level.

Fentianyl, belongs to powerful narcotics analgesics, and is a µ special opiates receptor stimulant. Fentanyl is one of the varieties that been listed in management of United Nations "Single Convention of narcotic drug in 1961". Among the opiates agents that under international control, fentanyl is one of the most commonly used to cure moderate to severe pain1. After continuous injection of fentanyl, the sufferer will have the performance of protracted opioid abstinence syndrome, such as ataxia and irritability etc2,3, which presents the addiction after taking fentanyl in a long time. Compared with drug addicts of amphetamine, drug addicts who take fentanyl mainly have got the possibility of higher infection rate of HIV, more dangerous injection behavior and more lifelong medication overdose 4.

more dangerous injection behavior and more lifelong medication overdose 4.

The Multi-Drug Rapid Test Panel is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of FYL in urine. The Multi-Drug Rapid Test Panel yields a positive result when FYL in urine exceeds detective level.

Propoxyphene (PPX) is a narcotic analgesic compound bearing structural similarity to methadone. As an analgesic Propoxyphene (2ra) is a narcoit enangesic compound bearing structural similarity to metinatione. As an analgesic, Propoxyphene can be from 50-75% as potent as oral codeine. Darvocet<sup>18</sup>, one of the most common brand names for the drug, contains 50-100 mg of Propoxyphene apsylate and 325-650 mg of acetaminophen. Peak plasma concentrations of Propoxyphene are achieved from 1 to 2 hours post dose. In the case of overdose, Propoxyphene blood concentrations can reach significantly higher levels. In humans, Propoxyphene is metabolized by N-demethylation to yield Norpropoxyphene has a longer half-life (30 to 36 hours) than parent Propoxyphene (6 to 12 hours). The accumulation of Norpropoxyphene seen with repeated doses may be largely responsible for resultant toxicity.

The Multi-Drug Rapid Test Panel yields a positive result when the concentration of Propoxyphene in urine exceeds

### Synthetic Cannabis(K2)

Synthetic Cannabis (T2)
Synthetic Cannabis or K2 a psychoactive herbal and chemical product that, when consumed, mimics the effects of Marijuana. It is best known by the brand names K2 and Spice, both of which have largely become genericized trademarks used to refer to any synthetic Cannabis product. The studies suggest that synthetic Cannabis intoxication is associated with acute psychosis, worsening of previously stable psychotic disorders, and also may have the ability to trigger actronic (long-term) psychotic disorder among vulnerable individuals such as those with a family history of mental illness. Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 72 hours after smoking (depending on usage/dosage). As of March 1, 2011, five cannabinoids, JWH-018, JWH-073, CP-47, JWH-200and cannabicyclobexanol are now illegal in the US because these substances have the potential to be extremely harmful and, herefore, pose an imminent hazard to the public safety. The Multi-Drug Rapid Test yields a positive result when the synthetic Cannabis metabolite in urine exceeds detective level.

Prevabaling/GBB

Pregabalin(PGB)
PGB regulates the influx of calcium ions in neurons and reduces the release of excitatory neurotransmitters such as glutamate, norepinephrine, and substance P, thereby inhibiting neuronal excitability, especially in the cerebral cortex, amygdala, and hippocampus. An area of the nervous system that is rich in synaptic connections. It is generally believe PGB is well tolerated, but in recent years, reports of PGB intolerance have increased year by year, which has aroused people's attention. Some cases suggest that PGB may have been abused, especially in substance-dependent individuals most of the current PGB abuse and addictions are caused by excessively increasing PGB doses or in combination with

drug
The Multi-Drug Rapid Test Panel yields a positive result when the concentration of Propoxyphene in urine exceeds

6-Monoacetylmorphine(6-MAM)

6-Monoacetylmorphine(6-MAM)

6-Monoacetylmorphine(6-MAM) or 6-acetylmorphine(6-AM)is one of three active metabolites of heroin(diacetylmorphine), the others being morphine and the much less active 3-monoacetylmorphine(3-MAM).6-MAM is rapidly created from heroin in the body, and then is either metabolized into morphine or excreted in the urine.6-MAM remains in the urine for no more than 24 hours. So a urine specimen must be collected soon after the last heroin use, but the presence of 6-MAM guarantees that heroin was in fact used as recently as within the last day, 6-MAM is naturally found in the brain, but in such small quantities that detection of this compound in urine virtually guarantees that heroin has recently been consumed. The Multi-Drug Rapid Test Panel yields a positive result when the concentration of 6-Monoacetylmorphine in urine exceeds detective level.

detective level.

Methcathinone(MCT)

Metheathinone, is a monoamine alkaloid and psychoactive stimulant, a substituted cathinone. Metheathinone is a highly addictive drug, primarily psychologically addicting and most of the signs of addiction to the drug are emotional or psychological. It has been popularized and continues to be sold under misleading names such as "bath salts", "plant fertilizers" or "research chemicals", but it is actually a powerful psycho-stimulant used as a recreational drug. Effects of this drug typically last from 4 to 6 hours. It is used as a recreational drug due to its potent stimulant and euphoric effects and is considered to be addictive, with both physical and psychological withdrawal occurring if its use is discontinued after prolonged or high-dosage administration replaced to be added to the country with own physical and psychological winding in a set of discontinuous prolonged or high-dosage administration. The Multi-Drug Rapid Test Panel yields a positive result when the concentration of Methcathinone in urine exc

## Ethyl- B-D-Glucuronide(ETG)

Ethyl-β-D-Glucuronide(ETG)
Ethyl Glucuronide (ETG) is a metabolite of ethyl alcohol which is formed in the body by glucuronidation following exposure to tehanol, such as by drinking alcoholic beverages. It is used as a biomarker to test for ethanol use and to mo alcohol abstinence in situations where drinking is prohibited, such as in the military, in professional monitoring programs(health professionals, attorneys, airline pilots inrecovery from addictions), in schools, in liver transplant clinic in recovering alcoholic patients. ETG can be measured in urine up to approximately 80 hours after ethanol is ingested. It is a more accurate indicator of the recent exposure to alcohol than measuring for the presence of ethanol itself. The Multi-Drug Rapid Test yields a positive result when the concentration of Ethyl Glucuronide inurine exceeds detective I IPRECAUTIONS!

- For healthcare professionals including professionals at point of care sites.

  Immunoassay for *In vitro* diagnostic use only. The test Panel should remain in the sealed pouch until use.
  All specimens should be considered potentially hazardous and handled in the same manner as an infectiou
  The used test Panel should be discarded according to federal, state and local regulations.
- s drug-bovine protein antigen conjugate on the membrane and the conjugate pad of each test contains

## [STORAGE AND STABILITY]

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test Panels must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

ISPECIMEN COLLECTION AND PREPARATION1 The urine specimen should be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed well before testing. When testing cards with S.V.T. or Alcohol storage of urine specimens should not exceed 2 hours at room temperature or 4 hours refrigerated prior to

Waterials Frovided								
Test Panel     Package insert								
Adulteration Color Chart (when applicable)								
Materials	Required But Not Provided							
Specimen collection container	• timer							

## Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

1.Bring the pouch to room temperature before opening it. Remove the test panel from thesealed pouch and use it as soon as

possible.

2. Take off the cap outside of the test end. With arrows pointing toward the urine specimen, immerse the test panel vertically into the urine specimenfor at least 8-10 seconds. Immerse the test panel to at least the level of the wavy lines on the strip(s) do not pass the arrows on the test panel when immersing the panel.

3.Place the test panel on a non-absorbent flat'surface, start the timer and wait for the red line(s) to appear. The results should be read at 5 minutes. Do not interpret results after 10 minutes.

## [INTERPRETATION OF RESULTS]

POSITIVE RESULT:	Only one colored band appears in the control region (C).  No apparent colored band appears in the test region (T).
NEGATIVE RESULT:	Two colored bands appear on the membrane. One band appears in the control region (C) and another band appears in the test region (T).
INVALID RESULT:	Control band fails to appear. Results from any test which has not produced a control band at the specified reading time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

A procedural control is included in the test. A line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LLIMITATIONS)

1. The Multi-Drug Rapid Test Panel provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory 2. There is a possibility that technical or procedural errors, as well as interfering substances in the urine specimen may

erroneous results.

3.Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.

4.A positive result does not indicate level or intoxication, administration route or concentration in urine.

5.A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.

below the cut-off level of the test.

6.This test does not distinguish between drugs of abuse and certain medications.

7.A positive test result may be obtained from certain foods or food supplements. Alcohol in the atmosphere, such as spray from perfumes, deodorizers, glass cleaners etc. can affect the Alcohol Rapid Tests. Therefore, adequate measures should taken to avoid undue interference from such atmospheric agents in the testing area.

8.The test is only for detection of presence/ absence of alcohol in the urine, which may result from habitual drinking or medications and does not discriminate the two.

TEXPECTED VALUES]

The negative result indicates that the drug concentration is below the detectable level. Positive result means the concentration of drug is above the detectable level.

## [PERFORMANCE CHARACTERISTICS]

Mothod

A side-by-side comparison was conducted using the Multi-Drug Rapid Test Panel and commercially available drug rapid tests. Testing was performed on approximately 250 specimens per drug type previously collected from subjects presenting for Drug Screen Testing. Presumptive positive results were confirmed by GC/MS.

CC/MS

Method			GC/MS	% agreement with GC/MS		
Multi-Drug R	Lapid Test Panel	Positive	Negative			
AMP 300	Positive	105	0	>99.9%		
AMP 300	Negative	0	145	>99.9%		
BAR 300	Positive	102	0	>99.9%		
DAK 300	Negative	0	148	>99.9%		
BZO 300	Positive	123	0	>99.9%		
BZO 300	Negative	0	127	>99.9%		
DUD 10	Positive	106	0	>99.9%		
BUP 10	Negative	0	144	>99.9%		
COC 150	Positive	113	0	>99.9%		
COC 150	Negative	0	137	>99.9%		
THC 50	Positive	94	0	>99.9%		
THC 30	Negative	0	156	>99.9%		
MTD 300	Positive	90	0	>99.9%		
M1D 300	Negative	0	160	>99.9%		
MET 300	Positive	79	0	>99.9%		
MET 300	Negative	0	171	>99.9%		
MDM 1 200	Positive	104	0	>99.9%		
MDMA 300	Negative	0	146	>99.9%		
MOP 300	Positive	100	0	>99.9%		
MOP 300	Negative	0	150	>99.9%		
OBI 2000	Positive	121	0	>99.9%		
OPI 2000	Negative	0	129	>99.9%		
non as	Positive	92	0	>99.9%		
PCP 25	Negative	0	158	>99.9%		
TCA 1000	Positive	96	0	>99.9%		
TCA 1000	Negative	0	154	>99.9%		

I	Method		GC/MS	% agreement with GC/MS
Multi-Drug I	Rapid Test Panel	Positive	Negative	
TRA 100	Positive	93	0	>99.9%
1 KA 100	Negative	0	157	>99.9%
KET 300	Positive	79	0	>99.9%
KE1 500	Negative	0	171	>99.9%
OXY 100	Positive	86	0	>99.9%
OX 1 100	Negative	0	164	>99.9%
EDDD 100	Positive	98	0	>99.9%
EDDP 100	Negative	0	152	>99.9%
FMT 200	Positive	80	0	>99.9%
FYL 200	Negative	0	170	>99.9%
DDM 200	Positive	96	0	>99.9%
PPX 300	Negative	0	154	>99.9%
152.50	Positive	79	0	>99.9%
K2 50	Negative	0	159	>99.9%
DCD 500	Positive	92	0	>99.9%
PGB 500	Negative	0	169	>99.9%
	Positive	99	0	>99.9%
6-MAM 10	Negative	0	157	>99.9%
	Positive	82	0	>99.9%
MCT 500	Negative	0	172	>99.9%
EMG 400	Positive	98	0	>99.9%
ETG 500	Negative	0	165	>99.9%

			% A	greeme	nt with	Comme	rcial Ki	t				
	AMP 300	BAR 300	BZO 300	BUP 10	COC 150	THC 50	MTD 300	MET 300	TRA 100	MDMA 300	MOP 300	OPI 2000
Positive Agreement	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%
Negative Agreement	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%
Total Results	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%
	PCP 25	TCA 1000	KET 300	OXY 100	EDDP 100	FYL 200	PPX 300	K2 50	PGB 500	6-MAM 10	MCT 500	ETG 5000
Positive Agreement	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%
Negative Agreement	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%
Total Results	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%	>99.9%

A study was conducted at three hospitals by Jaypersons using three different lots of product to demonstrate the within run between run and between operator precision. An identical card of coded specimens, containing drugs at concentrations of  $\pm$  50% and  $\pm$  25% cut-off level, was labeled, blinded and tested at each site. The results are given below:

Site A

Site B

## AMPHETAMINE (AMP 300)

	n non site	Sit	te A	Site	В	Site C		
Amphetamine conc. (ng/mL)	n per site	-	+	-	+	-	+	
0	10	10	0	10	0	10	0	
150	10	10	0	10	0	10	0	
225	10	9	1	8	2	9	1	
375	10	1	9	2	8	2	8	
450	10	0	10	0	10	0	10	
BARBITURATES (BAR 300)			•					
		Sit	te A	Site B		Site C		
Secobarbital conc. (ng/mL)	n per site	-	+	-	+	-	+	
0	10	10	0	10	0	10	0	
150	10	10	0	10	0	10	0	
225	10	9	1	8	2	9	1	
375	10	2	8	1	9	2	- 8	

	Site A		Site	В	Site C	
n per site	-	+	-	+	-	+
10	10	0	10	0	10	0
10	10	0	10	0	10	0
10	9	1	9	1	9	1
10	1	9	1	9	1	9
10	0	10	0	10	0	10
	10 10 10	n per site - 10 10 10 10 10 10 10 10 10 10 10 10 10	n per site	n per site	n per site - + - + 10 10 0 10 0 10 10 10 0 10 10 9 1 9 1 10 1 9 1 9	n per site         -         +         -         +         -           10         10         0         10         0         10           10         10         0         10         0         10           10         9         1         9         1         9         1           10         1         9         1         9         1         9         1

## Buprenorphine (BUP 10)

D 11 ( 1 T)			*				
Buprenorphine conc. (ng/mL)	n per site	-	+	-	+	-	+
0	10	10	0	10	0	10	0
5	10	10	0	10	0	10	0
7.5	10	9	1	9	1	8	2
12.5	10	1	9	1	9	1	9
15	10	0	10	0	10	0	10

Benzoylecgonine conc. (ng/mL)	1						
0	n per site	31	te A	Site	+	Site	C +
0	10	10	0	10	0	10	0
75	10	10	0	10	0	10	0
112.5	10	9	1	9	1	9	1
187.5	10	1	9	1	9	1	9
225 ARIJUANA (THC50)	10	0	10	0	10	0	10
		Si	te A	Site	В	Site	С
11-nor-△9-COOH conc. (ng/mL)	n per site		+		+		+
0	10	10	0	10	0	10	0
25	10	10 9	0	10	0	10 9	0
37.5 62.5	10	1	1 9	8	9	2	8
75	10	0	10	0	10	0	10
ETHADONE (MTD300)							
Methadone conc. (ng/mL)	n per site	Si	te A	Site		Site	
		4.0	+	10	+	10	+
150	10	10	0	10	0	10 10	0
225	10	9	1	9	1	9	1
375	10	1	9	1	9	1	9
450	10	0	10	0	10	0	10
ETHAMPHETAMINE (MET300)							
Methamphetamine conc. (ng/mL)	n per site	Si	te A	Site		Site	
0	10	10	+ 0	10	+	10	+ 0
150	10	10	0	10	0	10	0
225	10	9	1	9	1	9	1
375	10	1	9	2	8	1	9
450	10	0	10	0	10	0	10
ETHYLENEDIOXYMETHAMPHETA	MINE (MDM/	_			P		
Methylenedioxymethamphetamine conc. (ng/mL)	n per site	Si	te A	Site	B +	Site	+
0	10	10	0	10	0	10	0
150	10	10	0	10	0	10	0
225	10	8	2	9	1	9	1
375	10	1	9	1	9	1	9
450	10	0	10	0	10	0	10
ORPHINE (MOP 300)		Si	te A	Site	В	Site	С
Morphine conc. (ng/mL)	n per site		+		+		+
0	10	10	0	10	0	10	0
150	10	10	0	10	0	10	0
225	10	9	1	9	1	9	1
375 450	10	0	9	0	9	0	9 10
ORPHINE/OPIATE (OPI 2000)	10	0	10		10		10
	T .	Si	te A	Site	В	Site	С
Morphine conc. (ng/mL)	n per site		+		+		+
0	10	10	0	10	0	10	0
1000	10	10 9	0	10 9	0	10 9	0
1500 2500	10	1	9	1	9	1	1 9
3000	10	0	10	0	10	0	10
ENCYCLIDINE (PCP 25)							
Phencyclidine conc. (ng/mL)	n per site	Si	te A	Site		Site	
			+		+		+
0 12.5	10	10	0	10	0	10 10	0
12.5	10	8	2	9	1	9	1
31.25	10	1	9	1	9	1	9
37.5	10	0	10	0	10	0	10
ICYCLIC ANTIDEPRESSANTS (TC.	A 1,000)						
Nortriptyline conc. (ng/mL)	n per site	Si	te A	Site		Site	
0		10	+ 0	10	+ 0	10	+ 0
500	10	10	0	10	0	10	0
750	10	9	1	8	2	8	2
1,250	10	1	9	1	9	1	9

COCAINE (COC 150)

MARIJUANA (THC50)

METHADONE (MTD300)

METHAMPHETAMINE (MET300)

METHYLENEDIOXYMETHAMPHETAMINE (MDMA300) Ecstasy

MORPHINE (MOP 300)

MORPHINE/OPIATE (OPI 2,000)

PHENCYCLIDINE (PCP25)

TRICYCLIC ANTIDEPRESSANTS (TCA1000)

TRAMADOL (TRA 100)

KETAMINE (KET300)

Ecgonine

△%-THC

(±) -3,4-Methyle

Mephentermine

amphetamine

Oxycodone

Thebaine

Norcodeine

Normorphon

Oxymorphon

Oxycodone

Procaine

Imipramine Clomipramir

Doxepine

Promethazin

Perphenazine

Phencyclidine

EDDP

MDE

d,l-O-Desmethyl venlafaxine

(+) Chlorpheniramine

4-Hydroxyphencyclidine

d-Methamphetamin

l-Methamphetamin

Thioridazine

Dithiaden

methamphetamine

3,4-Methylenedioxyethyl

200

25,000

300

300

1,500

6,000

50,000

50,000

15,000

3,000

1,000

3,000

3,000

2,000

200

100.000

2,000

25,000

25,000

25,000

25,000

25,000

25,000

100,000

25,000

500

20,000

Benzoylecgonine Cocaine HCl

Methadon

11-nor-△8-THC-9 COOH

ρ-Hydroxymethamphetamine

D-Methamphetamine

L-Methamphetamine

(±) 3,4-Methylenediox

methamphetamine HCl

(±) 3,4-Methylenedioxy amphetamine HCl

Morphine-3-β-D-Glucuronide

Hydrocodone

Hydromorphor

Ethylmorphine Hydrocodone

Hydromorphone

Phencyclidine

Nortriptyline

Trimipramine

Amitriptyline

Desipramine

Cyclobenzaprine

Cis-tramadol

Procyclidine

Promazine

Promethazin

Pentazocine

Phencyclidine

Disopyramide

Tetrahydrozolin

(1R, 2S) - (-)-Ephedrine

Dextromethorphan

Methoxyphenamine

n-Desmethyl-cis-tramadol

6-Monoacetylmorphine

Morphine 3-β-D-glucuronio

Codeine

6-Monoacethylmorphine

n per site  10  10  10  10  10  10  10  10  10	0 Sit 10 10 9 1	10 te A + 0	0 Site	10	0	10	450 Oxycodone (OXY100)	10	0	10	0	10	0
10 10 10 10 10 10	10 10 9	+ 0	Site			_	Oxycodone (OXY100)						
10 10 10 10 10 10	10 10 9	+ 0	Site										
10 10 10 10 10 n per site	10 9 1	_		B +	Site	+ +	Oxycodone conc. (ng/mL)	n per site	Sit	te A	Site	B +	١.
10 10 10 n per site	9		10	0	10	0	0	10	10	0	10	0	10
10 10 n per site	1	0	10	0	10	0	50	10	10	0	10	-0	10
n per site		1	9	1	9	1	75	10	9	1	9	1	9
n per site	0	9	1	9	1	9	125	10	1	9	1	9	1
10		10	0	10	0	10	150	10	0	10	0	10	0
10							2-Ethylidene-1,5-dimethyl-3,3-diphenylpyr	rolidine (EDD					
	Sit	te A	Site	B +	Site	+ +	EDDP conc. (ng/mL)	n per site	Sit	te A	Site	B +	<del> </del>
	10	0	10	0	10	0	0	10	10	0	10	0	10
10	10	0	10	0	10	0	50	10	10	0	10	0	10
10	9	1	9	1	9	1	75	10	9	1	9	1	9
10	1	9	2	8	1	9	125	10	1	9	1	9	1
10	0	10	0	10	0	10	150	10	0	10	0	10	0
AINE (MDM/	1300) Ecs	tasy					Fentanyl (FYL200)						
n nor cito	Sit	e A	Site		Site	C	EVI cone (ng/mI)	n ner site	Sit				_
		+	<u> </u>	+		+		-	-		-		_
		_	_		_	_							10
10		0	10	0	10	0							10
10	8	2	9	1	9	1			9				9
10	1	9	1	9	1	9			1	_			1
10	0	10	0	10	0	10		10	0	10	0	10	0
							Propoxyphene (PPX300)						
n per site	Sit	te A	Site	B +	Site	+ +	PPX conc. (ng/mL)	n per site	Sit	te A	Site	B +	-
10	10	0	10	0	10	0	0	10	10	0	10	0	10
10	10	0	10	0	10	0	150	10	10	0	10	0	10
10	9	1	9	1	9	1	225	10	9	1	9	1	9
10	1	9	1	9	1	9	375	10	1	9	1	9	1
10	0	10	0	10	0	10	450	10	0	10	0	10	0
							Synthetic Cannabis(K2 50)						
	Sit	le A	Site	В	Site	e C			Sit	te A	Site	В	
n per site		+		+		+	Synthetic Cannabis conc. (ng/mL)	n per site	-	+	-	+	Τ.
10	10	0	10	0	10	0	0	10	10	0	10	0	10
10	10	0	10	0	10	0	25	10	10	0	10	0	10
10	9	1	9	1	9	1	37.5	10	9	1	9	1	9
10	1	9	1	9	1	9	62.5	10	1	9	2	8	1
10	0	10	0	10	0	10	75						
	$\overline{}$			•	•			10	0	10	0	10	0
							Pregabalin(PGB 500)	10	0	10	0	10	0
	Sit	te A	Site	В	Site	· C				te A	Site		0
n per site	Sit	te A	Site	B +	Site	+ +	Pregabalin(PGB 500)  Pregabalin conc. (ng/mL)	n per site					0
n per site	Sit		Site 10		Site				Sit	te A	Site	В	
		+		+		+	Pregabalin conc. (ng/mL)	n per site	Sit	te A	Site	B +	10
10	10	+ 0	10	+ 0	10	+ 0	Pregabalin conc. (ng/mL)	n per site	Sin	te A + 0	Site	B + 0	10
10	10 10	+ 0 0	10 10	+ 0	10	+ 0 0	Pregabalin conc. (ng/mL)  0  250	n per site 10 10	10 10	te A + 0 0	Site - 10 10	B + 0 0	10
10 10 10	10 10 8	+ 0 0 2	10 10 9	+ 0 0	10 10 9	+ 0 0	Pregabalin conc. (ng/mL)  0  250  375	n per site 10 10 10	10 10 9	te A	Site - 10 - 10 - 9	B + 0 0 1	10 10 10 9 1 1
10 10 10 10	10 10 8 1	+ 0 0 2 9	10 10 9	+ 0 0 1 9	10 10 9 1	+ 0 0 1 9	Pregabalin conc. (ng/mL)  0  250  375  625	n per site  10  10  10  10	10 10 10 9 1	te A + 0 0 1 1 9 10	Site	B + 0 0 1 1 8 10	10 10 9 1
10 10 10 10 10 10 1,000)	10 10 8 1	+ 0 0 2 9	10 10 9	+ 0 0 1 9 10	10 10 9	+ 0 0 1 9	Pregabalin conc. (ng/mL)  0 250 375 625 750 6-Monoacetylmorphine(6-MAM 10)	n per site  10  10  10  10  10	Sit 10 10 9 1 0 Sit	te A	Site - 10 10 9 2	B + 0 0 1 8 10 B	10 10 9 1
10 10 10 10 10 10 1,000)	10 10 8 1 0	+ 0 0 2 2 9 10 te A +	10 10 9 1 0	+ 0 0 1 9 10 B	10 10 9 1 0	+ 0 0 0 1 1 9 10 C C +	Pregabalin conc. (ng/mL)  0 250 375 625 750 6-Monoacetylmorphine(6-MAM 10)  6-MAM conc. (ng/mL)	n per site  10 10 10 10 10 10 n per site	Sin - 10 10 9 1 0 Sin	te A	Site   -   10   10   9   2   0     Site   -	B + 0 0 1 1 8 10 B + +	10 10 10 9 9 1 1 0
10 10 10 10 10 1,000)	10 10 8 1 0	+ 0 0 0 2 9 10 te A + 0	10 10 9 1 0 Site	+ 0 0 0 1 1 9 10 B + 0	10 10 9 1 0	+ 0 0 0 1 1 9 10 C C + 0	Pregabalin cone. (ng/mL)  0 250 375 625 750 6-Monoacetylmorphine(6-MAM 10)  6-MAM cone. (ng/mL)  0	n per site  10  10  10  10  10  10  n per site  10	Sit 10 10 9 1 0 Sit - 10	te A	Site - 10 Site - 10	B + 0 0 10 B + 0	10 10 10 9 1 1 0
10 10 10 10 10 10 1,000) n per site	10 10 8 1 0 Sit	+ 0 0 0 2 9 10 te A + 0 0 0	10 10 9 1 0 Site	+ 0 0 1 1 9 1 1 0 B + 0 0 0	10 10 9 1 0 Site	+ 0 0 0 1 1 9 10 C C + 0 0 0	Pregabalin conc. (ng/mL)  0 250 375 625 750 6-Monoacetylmorphine(6-MAM 10)  6-MAM conc. (ng/mL)  0 5	n per site  10  10  10  10  10  10  10  10  10  1	Sin 10 10 9 1 0 Sin - 10 10 10	te A	Site - 10 10 10 10 10 10 10 10 10 10 10 10 10	B + 0 0 1 1 8 10 B + 0 0 0	10 10 9 11 0
10 10 10 10 10 10 1,000) n per site 10 10	10 10 8 1 0 Sit	+ 0 0 0 2 9 10 te A + 0 0 1 1	10 10 9 1 0 Site	+ 0 0 1 9 10 B + 0 0 0 2	10 10 9 1 0 Site	+ 0 0 0 1 1 9 10 CC + 0 0 0 2	Pregabalin conc. (ng/mL)  0 250 375 625 750 6-Monoacetylmorphine(6-MAM 10)  6-MAM conc. (ng/mL)  0 5 7.5	n per site  10  10  10  10  10  10  10  10  10  1	Sit - 10 10 9 1 1 0 Sit - 10 10 10 9 9 1 10 10 9 9	te A	Site   -   10   10   9   2   0     Site   -   10   10   10   10   10   10   9   10   10	B + 0 0 1 1 B + 0 0 0 1 1	10 10 9 11 10 0
10 10 10 10 10 10 1,000) n per site	10 10 8 1 0 Sit	+ 0 0 0 2 9 10 te A + 0 0 0	10 10 9 1 0 Site	+ 0 0 1 1 9 1 1 0 B + 0 0 0	10 10 9 1 0 Site	+ 0 0 0 1 1 9 10 C C + 0 0 0	Pregabalin conc. (ng/mL)  0 250 375 625 750 6-Monoacetylmorphine(6-MAM 10)  6-MAM conc. (ng/mL)  0 5	n per site  10  10  10  10  10  10  10  10  10  1	Sin 10 10 9 1 0 Sin - 10 10 10	te A	Site - 10 10 10 10 10 10 10 10 10 10 10 10 10	B + 0 0 1 1 8 10 B + 0 0 0	10 10 9 11 0
	n per site  10  10  10  10  10  10  10  10  10  1	Site	10	Site A   Site	Site A   Site B	Site A   Site B   Site	Site A   Site B   Site C	n per site         Site A         Site B         Site C         FYL conc. (ng/mL)           10         10         0         10         0         10         0           10         10         0         10         0         100         0         100           10         8         2         9         1         9         1         9         1         100         100         100         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         300         100         100         300         100         100         100         100         1150         00         1150         1150         1150         1150         1150         1150         1150	n per site         Site A         Site B         Site C         FYL conc. (ng/mL)         n per site           10         10         0         10         0         10         0         10         0         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10 </td <td>  N per site   Site A</td> <td>  N per site   Site  </td> <td>  New Site   S</td> <td>n per site         Site A         Site B         Site C         FYL conc. (ng/mL)         n per site         C         A         Site B         Site B         Site B         PPL conc. (ng/mL)         n per site         C         A         C         A         C         A         A         C         A         A         C         A         A         C         A         A         C         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         B         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D</td>	N per site   Site A	N per site   Site	New Site   S	n per site         Site A         Site B         Site C         FYL conc. (ng/mL)         n per site         C         A         Site B         Site B         Site B         PPL conc. (ng/mL)         n per site         C         A         C         A         C         A         A         C         A         A         C         A         A         C         A         A         C         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         A         B         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D         D

30,000

17.000

17,000

100,000

50,000

300

50,000

30,000

50,000

15,000

300

50,000

25,000 25,000

50,000 25,000

12,500

50.000

2,000

2,000

50,000

10,000

10,000

50,000

25,000

25,000

50,000

50,000

50,000

50,000

50,000

50,000

100,000

50,000

100,000

TRAMADOL (TRA 100)

Note		p	-	+	-	+	-	+
So	0	10	10	0	10	0	10	0
								0
125				_	_			
No.   No.					_			2
	125	10	1	9	1	9	2	8
Negation   Negation	150	10	0	10	0	10	0	10
Negation   Negation	KETAMINE (KET300)	•	•	•				
Netamine cone. (ng/mL)	·		Sit	e A	Site	В	Site	C
	Ketamine conc. (ng/mL)	n per site	_					+
150								
			10	0	10	0	10	0
375   10	150	10	10	0	10	0	10	0
375   10	225	10	9	1	8	2	9	1
			1				2	8
Daysedone (DXY100)				_				
Site A   Site B   Site C		10	0	10	0	10	0	10
Oxycodone conc. (ng/mL)	Oxycodone (OXY100)							
			Sit	te A	Site	В	Site	C
0	Oxycodone conc. (ng/mL)	n per site	-	+	-	+	-	+
Simple	0	10	10		10	0	10	0
75								0
125								
150	75		9			1	9	1
Pethylidene-1,5-dimethyl-3,3-diphenylpyroldine (EDDP)   EDDP conc. (ng/mL)	125	10	1	9	1	9	1	9
Pethylidene-1,5-dimethyl-3,3-diphenylpyroldine (EDDP)   EDDP conc. (ng/mL)	150	10	0	10	0	10	0	10
Bernome			P100)					
EDDP conc. (ng/mL)	- Etnyhaene 1,5 annetnyr 5,5 aipnenyrpy	Tronume (EDD		. A	Cita	D	Cita	C .
0	EDDP conc. (ng/mL)	n per site						
Simple						+		+
75	0	10	10	0	10	0	10	0
75	50	10	10	0	10	0	10	0
125	75			1		1		1
The state of the							_	9
FYL conc. (ng/mL)								
FYL conc. (ng/mL)		10	- 0	10	0	10	0	10
PYL conc. (ng/mL)	Fentanyl (FYL200)							
10			Sit	te A	Site	В	Site	C
100	FYL conc. (ng/mL)	n per site	-	+	-	+	-	+
100	0	10	10	0	10	0	10	0
150								
250								0
PPX conc. (ng/mL)	150	10	9	1	9	1	9	1
Propagation	250	10	1	9	1	9	1	9
Propagation	300	10	0	10	0	10	0	10
PPX conc. (ng/mL)								
PPX conc. (ng/mL)	roponypiene (rriness)		e:	10 A	Cito	D	Cita	C
10	PPX conc. (ng/mL)	n per site						_
150								+
225	0	10	10	0	10	0	10	0
375	150	10	10	0	10	0	10	0
375	225	10	9	1	9	1	9	1
Synthetic Cannabis (KZ 50)								9
Synthetic Cannabis (Cannabis conc. (ng/mL)								
Synthetic Cannabis conc. (ng/mL)		10	- 0	10	0	10	0	10
Synthetic Cannabis cone. (ng/mL)	Synthetic Cannabis(K2 50)							
10			Sit	te A	Site	В	Site	C
0	Synthetic Cannabis conc. (ng/mL)	n per site	-	+	-	+	-	+
25	0	10			10		10	0
37.5   10   9   1   9   1   9   1   9   1   9   1   62.5   10   1   9   2   8   1   1   75   10   0   10   0   10   0   10   0   10   0								0
Company   Com								
75 10 0 10 0 10 0 10 0 10 0   Pregabalin (PGB 500)  Pregabalin conc. (ng/mL)								1
Pregabalin(PGB 500)  Pregabalin conc. (ng/mL)								9
Pregabalin conc. (ng/mL)		10	0	10	0	10	0	10
Pregabalin conc. (ng/mL)	Pregabalin(PGB 500)							
Pregabalin conc. (ng/mL)			Sit	te A	Site	В	Site	С
0	Pregabalin conc. (ng/mL)	n per site						+
250	0	10						0
375								
$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$								0
750	375	10	9	1	9	1	9	1
750	625	10	1	9	2	8	1	9
			0	10			0	10
6-MAM conc. (ng/mL)         n per site         Site A         Site B         Site C           0         10         10         0         10         0         10         0         10         0         10         0         10         0         10         0         10         0         10         10         0         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         10         <								10
6-MAM conc. (ng/mL)			en.	la A	eu-	D	014-	C
0 10 10 0 10 0 10 5 10 10 0 10 0 10 7.5 10 9 1 9 1 9 12.5 10 1 9 2 8 1	6-MAM conc. (ng/mL)	n per site						
5         10         10         0         10         0         10           7.5         10         9         1         9         1         9           12.5         10         1         9         2         8         1		-						+
7.5 10 9 1 9 1 9 12.5 10 1 9 2 8 1	0	10	10	0	10	0	10	0
7.5 10 9 1 9 1 9 12.5 10 1 9 2 8 1	5	10	10	0	10	0	10	0
12.5 10 1 9 2 8 1								1
								9
15   10   0   10   0   10   0								
	15	10	0	10	0	10	0	10

	Oxycodon	e (OXY100)	
Oxycodone	100	Hydromorphone	50,000
Oxymorphone	300	Naloxone	25,000
Levorphanol	50,000	Naltrexone	25,000
Hydrocodone	25,000		
2-Ethyliden	e-1,5-dimethyl-3,3-	diphenylpyrrolidine (EDDP100)	
2-Ethylidene-1,5-dimethyl-3,3-dipheny	lpyrrolidine (EDDI	?)	100
	Fentanyl	(FYL200)	
Alfentanyl	600,000	Buspirone	15,000
Fenfluramine	50,000	Fentanyl	100
Norfentanyl	20	Sufentanyl	50,000
Diazepam	300	Triazolam	5,000
Estazolam	1,250		
	Propoxyphe	ene (PPX300)	•
D-Propoxyphene	300	D-Norpropoxyphene	1500
	Synthetic Ca	nnabis(K2 50)	•
JWH-018 5-Pentanoic acid	50	JWH-018 5-Hydroxypentyl	500
JWH-073 4-butanoic acid	50	JWH-073 4-Hydroxybuty	500
JWH-018 4-Hydroxypentyl	400		
	Pregabali	n(PGB 500)	
Vigabatrin	>100000	Phenibut	100000
Gabapentin	>100000	Pregabalin	500
	6-Monoacetylmor	rphine(6-MAM 10)	•
Morphine	100000	6-Monoacetylmorphine	10
	Methcathino	one(MCT 500)	
S(-)-Methcathinone HCl	500	R(+)-Methcathinone HCl	1500
Methoxyphenamine	100000	3-Fluoromethcathinone HCl	1500
	Ethyl-β-D-Glucu	ronide(ETG 500)	•
Propyl β-D-glucuronide	50000	Glucuronic Acid	100000
Morphine 3β-glucuronide	100000	Ethanol	120000
Morphine 6β-glucuronide	100000	Methanol	130000
Ethyl- β-D-Glucuronide	500		

Effect of Urinary Specific Gravity
Fifteen (15) urine samples of normal, high, and low specific gravity ranges (1.005-1.045) were spiked with drugs at 50% below
and 50% above cut-off levels respectively. The Multi-Drug Rapid Test Panel was tested in duplicate using fifteen drug-free urine

and 50% above cut-off levels respectively. The Multi-Drug Rapid Test Panel was tested in duplicate using fifteen drug-free urine and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with drugs at 50% below and 50% above cut-off levels. The spiked, pH-adjusted urine was tested with the Multi-Drug Rapid Test Panel. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or drug positive urine containing, Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Marijunan, Methadone, Methadone, Methamphetamine, Morphine, Tramadol, Ketamine, Pheneyeldine, Tricyclic Antidepressants,Oxycodone, EDDP, Fentanyl, 3,4-Methylenedioxyamphetamine, Carfentanyl, Diazepam 3,4-methylenedioxyyamphetamine and Enthyl-4D-Glucuronide. The following compounds show no cross-reactivity when tested with the Multi-Drug Rapid Test Panel at a concentration of 100 µg/mL.

Non	Cross-Reacting	Compounds
Non	Cross-Reacting	Compound

Acetophenetidin	Cortisone	Zomepirac	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycorticosterone	Labetalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diflunisal	Isoxsuprine	Sulfamethazine
l-Ascorbic acid	Digoxin	d,l-Propanolol	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone
Atropine	β-Estradiol	Niacinamide	3-acetate
Benzilic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrocortisone
Benzoic acid	Erythromycin	Norethindrone	Tetrahydrozoline
Bilirubin	Fenoprofen	Noscapine	Thiamine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	Thioridazine
Caffeine	Gentisic acid	Oxalic acid	d,l-Tyrosine
Cannabidiol	Hemoglobin	Oxolinic acid	Tolbutamide
Chloral hydrate	Hydralazine	Oxymetazoline	Triamterene
Chloramphenicol	Hydrochlorothiazide	Papaverine	Trifluoperazine
Chlorothiazide	Hydrocortisone	Penicillin-G	Trimethoprim
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	d,l-Tryptophan
Chlorpromazine	3-Hydroxytyramine	Phenelzine	Uric acid
Cholesterol	d,l-Isoproterenol	Prednisone	Verapamil
Clonidine			

### Methcathinone(MCT 500)

Site 0

		Site A			В	Site C	
Methcathinone conc. (ng/mL)	n per site	-	+	-	+	-	+
0	10	10	0	10	0	10	0
250	10	10	0	10	0	10	0
375	10	9	1	9	1	9	1
625	10	1	9	2	8	1	9
750	10	0	10	0	10	0	10

### Ethyl- β-D-Glucuronide(ETG 500)

E4 1 0 B 61 11 ( 1 1 )		Sit	e A	Site	В	Site	C	
Ethyl- β-D-Glucuronide conc. (ng/mL)	n per site	-	+	-	+	-	+	
0	10	10	0	10	0	10	0	
250	10	10	0	10	0	10	0	
375	10	9	1	9	1	9	1	
625	10	1	9	2	8	1	9	
750	10	0	10	0	10	0	10	1

Analytical Sensitivity
A drug-free urine pool was spiked with drugs at the listed concentrations. The results are summarized below

Drug Concentration	AMI	300	BAR	300	BZO	300	BUP	10	COC	150	THC	50	MET	300	MDM	IA 300
Cut-off Range	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	26	4	27	3	27	3	26	4	26	4	26	4	27	3	25	5
Cut-off	15	15	16	14	15	15	14	16	13	17	14	16	16	14	14	16
+25% Cut-off	3	27	4	26	3	27	3	27	3	27	3	27	3	27	4	26
+50% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30
+300% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30
Drug Concentration	MOI	300	OPI	2000	PC	P 25	TCA	1000	OXY	100	MTI	300	TRA	100	KE	Т 300
Cut-off Range	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	27	3	27	3	25	5	25	5	27	3	26	4	27	3	27	3
Cut-off	15	15	14	16	15	15	15	15	15	15	14	16	15	15	15	15
+25% Cut-off	5	25	4	26	3	27	4	26	4	26	3	27	4	26	3	27
+50% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30
+300% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30
Drug Concentration	EDD	P 100	FYL	200	PPX	300	K2	50	PGE	3 500	6-MA	M 10	MCT	Г 500	ETG	5000
Cut-off Range	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	26	4	27	3	27	3	25	5	27	3	27	3	25	5	25	5
Cut-off	15	15	14	16	15	15	14	16	15	15	14	16	15	15	15	15
+25% Cut-off	3	27	4	26	3	27	4	26	5	25	4	26	3	27	4	26
+50% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30
+300% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30

Analytical Specificity
The following table lists the concentrations of compounds (ng/mL) that are detected as positive in urine by the Multi-Drug

Analytes	Concentration (ng/mL)	Analytes	Concentration (ng/mL)
	AMPHETAMIN	NE (AMP 300)	•
D,L-Amphetamine sulfate	1000	Phentermine	1,000
L-Amphetamine	25,000	Maprotiline	50,000
(±) 3,4-Methylenedioxy	500	Methoxyphenamine	6,000
amphetamine		D-Amphetamine	300
	BARBITURAT	ES (BAR 300)	
Amobarbital	5,000	Alphenol	600
5,5-Diphenylhydantoin	8,000	Aprobarbital	500
Allobarbital	600	Butabarbital	200
Barbital	8,000	Butalbital	8,000
Talbutal	200	Butethal	500
Cyclopentobarbital	30,000	Phenobarbital	300
Pentobarbital	8,000	Secobarbital	300
	BENZODIAZEP	INES (BZO300)	•
Alprazolam	100	Bromazepam	900
a-hydroxyalprazolam	1,500	Chlordiazepoxide	900
Clobazam	200	Nitrazepam	200
Clonazepam	500	Norchlordiazepoxide	100
Clorazepatedipotassium	500	Nordiazepam	900
Delorazepam	900	Oxazepam	300

[BIBLIOGRAPHY]

1. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

2. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735.

3. Stewart DJ, Inaba T, Lucassen M, Kalow W. Clin. Pharmacol. Ther. April 1979; 25 ed: 464, 264-8.

4. Ambre J. J. Anal. Toxicol. 1985; 9:241.

4.Ambre J. J. Anal. 10xtcol.1963; Y.241.

S. Winger, Gail, A Handbook of Drug and Alcohol Abuse, Third Edition, Oxford Press, 1992, page 146.

6.Robert DeCresce. Drug Testing in the workplace, 1989 page 114.

7.Glass, IB. The International Handbook of Addiction Behavior. Routledge Publishing, New York, NY. 1991; 216

7.Glass, IB. The International Handbook of Addiction Behavior, Routledge Publishing, New York, NY. 1991; 216

8.B. Cody, J.T., "Specimen Adulteration in drug urinalysis. Forensic Sci. Rev., 1990, 2:63.

9.C. Tsai, S.C. et.al., J. Anal. Toxicol. 1998; 22 (6): 474

10. Baselr R.C. Disposition of Toxic Drugs and Chemicals in Man. 6th Ed. Biomedical Publ., Foster City, CA 2002.

11. Hardman J.G., Limbird LE. Goodman and Gilman's: The Pharmacological Basis for Therapeutics. 10th Edition. McGraw Hill Medical Publishing, 2001; 208-209.

12. Volpicellim, Joseph R., M.D., Ph.D.: Alcohol Dependence: Diagnosis, Clinical Aspects and Biopsychosocial Causes, Substance Abuse Library, University of Pennsylvania, 1997.

33. \*Assexment of Zopiclone" (PDF). World Health Organization. Essential Medicines and Health Products World Health Organization, p. 9 (Section 5. Pharmacokinetics). Retrieved 5 December 2015.

14.Kratzsch C, Tenberken O, Peters FT et al. Screening, library-assisted identification, and validated quantification of 23

14.Kratzsch C, Jenberken O, Peters F I et al. Screening, library-assisted identification, and validated quantification of 23 benzodiazepines, funnazenii, saleplone, zoljeldem, and zojelonie in plasma by ljquid chromatography/mass spectrometry with atmospheric pressure chemical ionization. J. Mass Spec. 39: 856-872, 2004.
11. Sciatavsen I, Al-Sammurraie M, Mørland J, Brammes JG. Impairment related to blood drug concentrations of zopiclone and zolpidem compared with alcohol in apprehended drivers. Accid. Anal. Prev. 41: 462-466, 2009.
16.R. Baselt, Disposition of Toxic Drugs and Chemicals i Man, 8th edition, Biomedical Publications, Foster City, CA, 2008, pp. 1672-1679.

1677-1679.

17. Calkins RF, Aktan GB, Hussain KL (1995). "Methcathinone: the next illicit stimulant epidemic?". Journal of Psychoactive Drugs. 27 (3): 277-85. doi:10.1080/02791072.1995.10472472. PMID 8594170.

18. Methcathinone, https://en.wikipedia.org/wiki/Methcathinone.

19. Al-Motarreh, Ahmed, Baker, Kathryn; Broadley, Kenneth J. (2002). "Khat: Pharmacological and Medical Aspects and Its Social Use in Vennen". Phytotherapy Research 16 (2): 403-413. doi:10.1002/ptr.1106. PMID 12203257. Retrieved 11 March 2015 20.List of psychotropic substances under international control. International Narcotics Control Board. United Nations. Archived from the circuit on 2012.03.61

from the original on 2012-08-31.

from the original on 2012-08-31.

21.Hoffman, R; Al'Absi, M (December 2010). "Khat use and neurobehavioral functions: suggestions for future studies," (PDF). Journal of Ethiopharmacology 132 (3): 554-63. doi:10.1016/j.jep.2010.05.033. PMC 2976806. PMID 20555832

22."List of psychotropic substances under international control" (PDF). International Narcotics Control Board. Archived from the original (PDF) on 2012-08-31.

23.Bersani, F. S.; Corazza, O.; Simonato, P.; Mylokosta, A.; Levari, E.; Lovaste, R.; Schifano, F. (2013). "Drops of madness? Recreational misuse of tropicamide collyrium; early warning alerts from Russia and Italy". General Hospital Psychiatry 35 (5):571-3.Basel RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488 24.Malenka RC, Nestler EJ, Hyman SE (2009). "Chapter 15: Reinforcement and Addictive Disorders". In Sydor A, Brown RY. Molecular Neuropharmacology: A Foundation for Clinical Neuroscience (2nd ed.). NewYork: McGraw-Hill Medical. p. 375. ISBN 9780071481274.

25.American Psychiatric Association (2013). "Substance-Related and Addictive Disorders". American Psychiatric Publishing. pp 26. Juliano LM, Griffiths RR (2004). "A critical review of caffeine withdrawal: empirical validation of symptoms and signs

incidence, severity, and associated features". Psychopharmacology (Berl.) 176 (1):1–29. doi:10.1007/s00213-004-2000-x. PMID 15448977. Archived from the original on 29 January 2012.

27. Arnaud MJ. Pharmacokinetics and metabolism of natural methylxanthines in animal and man. Handb Exp Pharmacol 201

200:33-91.
28.Jeukendrup AE, Randell R-Fat burners: nutrition supplements that increase fat metabolism. Obes Rev 2011; 193:1-24.
29.Cumming, E. (22 April 2010). "Mephedrone: Chemistry lessons". London: The Daily Telegraph. Retrieved 2010-09-14.
30. "Drugs crackdown hailed a success". BBC News. 8 March 2010. Retrieved 2010-03-31.
31. Kihara, Rhiannon; Day, Edward (May 2014). "Transient psychotic episodes following recreational use of NRG-3". Progress in Neurology and Psychiatry 18 (3): 14–18. doi:10.1002/pnp.331. Retrieved 22 March2015.
32. Schiffano, F.; Albanese, A.; Fergus, S.; Sistin, J. L.; Deluca, P.; Cornazza, O.; Davey, Z.; Corkery, J.; Siemann, H.; Scherbaum, N.; Farre', M.; Torrens, M.; Demetrovics, Z.; Ghodse, A. H.; Psychonaut Web, M.; Rednet Research, G. (2010). "Mephedrone (4-methylmetheathinone; 'meow meow'): chemical, pharmacological and clinical issues". Psychopharmacology 214 (3):593–602. doi:10.1007/802013-010-2070-XISSN 0033-318S. PMID 21072502.
33. Work Group on Panic Disorder (January 2009). APA Practice Guideline for the Treatment of Patients With Panic Disorder (2nd ed.).

(2nd ed.).
34."FDA approved labeling for Xanax revision 08/23/2011" (PDF). Federal Drug Administration. 2011-08-23. p. 4. Retrieved

2011-09-14.
35."Xanax XR (Alprazolam) Clinical Pharmacology – Prescription Drugs and Medications". RxList. First DataBank, July 2008.

## Index of Symbols





Sejoy Biomedical Co., Ltd. Area C, Building 2, No.365, Wuzhou Road, Yuhang Economic Development Zone, Hangzhou City 311100 Zhejiang China





size: 300\*420mm